

Rehbar AI

Setup and User Guide

1. System Setup

Follow the steps below to install and prepare Rehbar AI on your laptop.

Step 1: Install Python

1. Download and install the latest version of Python.

Step 2: Open the Project Folder

1. Locate the folder containing the Rehbar AI project files.
2. Open the folder.

Step 3: Open Command Prompt in the Folder

1. Click on the address bar at the top of the folder window.
2. Type `cmd` and press **Enter**.
A Command Prompt window will open in that folder location.

Step 4: Install Required Libraries

In the Command Prompt window, copy and paste the following command, then press **Enter**:

```
pip install ultralytics opencv-python customtkinter pillow pyttsx3
```

Wait until all required packages are downloaded and installed. Once the installation completes without errors, the setup process is finished.

2. How to Use Rehbar AI

Step 1: Launch the Application

In the same Command Prompt window, type:

```
python main_rehbar.py
```

Press **Enter** to start the application.

Step 2: Start the System

1. Wait for the Rehbar dashboard to load.
2. Click the green **“START SYSTEM”** button.
3. The camera will activate, and you will hear a startup voice notification.

Step 3: Using the Camera

- Position the camera toward the area where navigation assistance is needed.
- The system will detect objects and people in view.
- When an object or person is recognized:
 - A bounding box will appear around them on the screen.
 - The system will announce their name through voice feedback.

To reduce unnecessary repetition, the system pauses briefly before repeating the name of the same object.

Step 4: Audio Settings

- Ensure your laptop speakers are turned on and set to an appropriate volume.
- Confirm that the **“Voice Feedback”** option is enabled on the dashboard.

Step 5: Stop the System

When you are finished:

1. Click the red **“STOP SYSTEM”** button.
2. The camera will turn off, and the application will close safely.

Your Rehbar AI system is now ready for use.